

Mapping of RNA modifications by direct Nanopore sequencing and JACUSA2

Christoph Dieterich^{*1,2,3}, Amina Lemsara^{1,2}, and Isabel Naarmann-de Vries^{1,2,3}

¹Klaus Tschira Institute for Integrative Computational Cardiology, University Heidelberg, 69120 Heidelberg, Germany

²Department of Internal Medicine III (Cardiology, Angiology, and Pneumology), University Hospital Heidelberg, 69120 Heidelberg, Germany

³German Centre for Cardiovascular Research (DZHK)-Partner Site Heidelberg/Mannheim, 69120 Heidelberg, Germany

Abstract

to be written

Keywords: Bayesian, 10X Genomics, Cell barcode assignment, Nonsense-mediated mRNA decay (NMD)

INTRODUCTION

Chemical modifications on DNA and histones, also known as epigenetics marks, strongly impact gene expression during cell differentiation and in several other biological programs. In the 1970s, it was recognized that RNA is also subjected to extensive covalent modification, and studies in the late 1980s revealed the widespread deamination of bases (termed RNA editing), which can lead to recoding if it occurs within coding sequences. Impressive development in the RNA modification field occurred during the past eight years, with the discovery of an extensive layer of base modifications in mRNAs. These can influence gene expression and have been already shown to be involved in primary cellular programs such as stem cell differentiation, response to stress, and the circadian clock. The study of RNA modifications and their effects is now referred to as epitranscriptomics, and it reveals striking similarities to what is known for epigenomics. To date thirteen distinct modifications have been identified on mRNA transcripts [?]. These modifications are catalyzed by a variety of dedicated enzymes and can be divided into two classes: modifications of cap-adjacent nucleotides and internal modifications.

*christoph.dieterich@uni-heidelberg.de

32 In contrast to the m7G cap, the impact of internal modifications on gene
 33 regulation has been less studied apart from RNA editing, which is mediated
 34 by RNA deaminases (e.g. the ADAR family). The most widespread internal
 35 mRNA modification is N6-methyladenosine (m6A). By modulating the pro-
 36 cessing of mRNA, m6A can regulate a wide range of physiological processes
 37 and its alteration has been linked to several diseases [???]. The modifica-
 38 tion is catalyzed co-transcriptionally by a Mega-Dalton methyltransferase
 39 complex, which includes the heterodimer METTL3-METTL14 and other
 40 associated subunits [?]. This modification is reversible since two proteins of
 41 the AlkB-family demethylases can remove m6A from mRNA transcripts [??].
 42 In mammals, m6A preferentially localizes within long internal exons and at
 43 the beginning of terminal exons at so-called DRACH motif (D = A/G/U,
 44 R = A/G, H = A/C/U) sites [???]. Once deposited, m6A is recognized
 45 by several reader proteins that can affect the fate of mRNA transcripts in
 46 nearly every step of the mRNA life cycle, which includes alternative splicing
 47 [?]. The best-described readers are the YTH domain family of proteins that
 48 decode the signal and mediate m6A functions. By affecting RNA structure,
 49 m6A can also indirectly influence the association of additional RNA-binding
 50 proteins (RBPs) and the assembly of larger messenger ribonucleoprotein
 51 particles (mRNPs).

52 Several approaches have been presented to map RNA modifications on
 53 RNA. Herein, we focus on mRNA modification site detection in general and
 54 on m6A in particular where antibody-based protocols (miCLIP), methylation-
 55 sensitive restriction enzyme assays (MazF) or transgenic approaches (TRIBE,
 56 DART) have been presented. All of the aforementioned approaches rely on
 57 high-throughput sequencing on the Illumina platform. This typically in-
 58 volves cDNA synthesis by reverse transcription and PCR-based library am-
 59 plification. One recent addition to the tool is direct RNA single molecule
 60 sequencing on the Oxford Nanopore Technology platform. While our software
 61 workflow is able to deal with Illumina and Nanopore-based approaches, the
 62 latter is the principal topic of our methods article.

63 MATERIALS

64 ONT direct RNA sequencing

- 65 1. 500 ng polyA⁺ RNA isolated from total RNA e.g. with Oligotex
 66 mRNA kit (Qiagen) or Dynabeads oligo dT₂₅ beads (Thermo Fisher
 67 Scientific) or *in vitro* transcriptome sample. Store RNA at -80 °C and
 68 the mRNA purification kit as recommended by the manufacturer.
- 69 2. Nuclease-free water. Store at room temperature.
- 70 3. Direct RNA-sequencing kit (SQK-RNA002, Oxford Nanopore Tech-
 71 nologies). Store at -20 °C.

- 72 4. NEBNext Quick Ligation Reaction Buffer (New England Biolabs).
73 Store at -20 °C.
- 74 5. T4 DNA Ligase (New England Biolabs). Store at -20 °C.
- 75 6. dNTP Mix (10 mM each). Store at -20 °C.
- 76 7. SuperScript IV Reverse Transcriptase (Thermo Fisher Scientific). Store
77 at -20 °C.
- 78 8. Agencourt RNAClean XP beads (Beckman Coulter). Store at 4 °C.
- 79 9. 70 % ethanol, freshly prepared.
- 80 10. Qubit dsDNA HS assay kit and Qubit Fluorometer (Thermo Fisher
81 Scientific).
- 82 11. Flow cell priming kit (EXP-FLP002, Oxford Nanopore Technologies).
83 Store at -20 °C.
- 84 12. Thermocycler.
- 85 13. Gentle rotator mixer.
- 86 14. Magnetic stand for 1.5 ml tubes.
- 87 15. 1.5 ml DNA LoBind tubes (Eppendorf), 0.2 ml PCR tubes.
- 88 16. MinION or GridION sequencing device and MinION R9.4.1 Flow cells
89 (FLO-MIN106D, Oxford Nanopore Technologies). Store Flow cells at
90 4 °C.

91 **Preparation of an *in vitro* transcriptome sample**

- 92 1. 100 ng polyA⁺ RNA isolated from total RNA e.g. with Oligotex
93 mRNA kit (Qiagen) or Dynabeads oligo dT₂₅ beads (Thermo Fisher
94 Scientific). Store RNA at -80 °C and the mRNA purification kit as
95 recommended by the manufacturer
- 96 2. 10 μM oligo(dT)-VN RT primer. TTTTTTTTTTTTTTTTTTTTTTTTTTTTTTVN.
97 Store at -20 °C.
- 98 3. 20 μM template switching oligo (TSO). ACTCTAATACGACTCAC-
99 TATAGGGAGAGGGCrGrG+G. Store at -20 °C.
- 100 4. 10 μM T7 extension primer. GCTCTAATACGACTCACTATAGG.
101 Store at -20 °C.
- 102 5. Nuclease-free water. Store at room temperature.

- 103 6. dNTP Mix (10 mM each). Store at -20 °C.
- 104 7. Template Switching RT Enzyme Mix (New England Biolabs). Store
105 at -20 °C.
- 106 8. Q5 Hot Start High-Fidelity 2X Master Mix (New England Biolabs).
107 Store at -20 °C.
- 108 9. RNase H (5,000 U/ml) (New England Biolabs). Store at -20 °C.
- 109 10. NucleoSpin Gel and PCR Clean-up, Mini kit for gel extraction and
110 PCR clean up (Macherey-Nagel) or equivalent. Store at room temper-
111 ature.
- 112 11. MEGAscript T7 transcription kit (Thermo Fisher Scientific). Store at
113 -20 °C.
- 114 12. RNA Clean & Concentrator-25 kit (Zymo Research). Store at room
115 temperature.
- 116 13. Thermocycler.
- 117 14. Table top centrifuge for 1.5 ml tubes.
- 118 15. Nanodrop spectrophotometer or equivalent.
- 119 16. 0.2 ml PCR tubes, 1.5 ml DNA LoBind tubes (Eppendorf).

120 **Hardware requirements**

121 All analyses have been performed/tested on two alternative hardware sys-
122 tems: a standard Linux desktop computer or an Apple iMac (Retina 5K,
123 ultimo 2014). The workflow requires a multi-core processor system with
124 minimal main memory of 16GB RAM and several GBs of free disk space
125 (depending on data set size).

126 **Software dependencies and installation**

127 Our analysis workflow has few requirements, which are detailed in Table 2.
128 Specifically, to execute our workflow, the following prerequisites are neces-
129 sary: a BASH shell, a JAVA runtime environment, a working PERL and
130 R installation. Additional i.e. non-standard software to process and map
131 Nanopore reads (bedtools, samtools and Minimap2) are obligatory, while
132 the installation of a Nanopore read simulator (NanoSim) is optional and de-
133 pends on your use case. Table ?? lists some additional R packages, which are
134 required to run the R code. Detailed instructions on how to setup are found
135 under https://github.com/dieterich-lab/MiMB_JACUSA2_chapter

136 METHODS

137 Overview Figure 1

138 Nanopore direct RNA sequencing

- 139 1. Adjust 500 ng polyA⁺ RNA to a total volume of 9 μ l with nuclease-
140 free water. Complete RT adapter ligation reaction (in 0.2 ml PCR
141 tube) with 3 μ l NEBNext Quick Ligation Reaction Buffer, 0.5 μ l
142 RNA CS (RCS, from SQK-RNA002), 1 μ l RT-Adapter (RTA, from
143 SQK-RNA002) and 1.5 μ l T4 DNA Ligase. Incubate 10 min at room
144 temperature.
- 145 2. Prepare reverse transcription master mix on ice during ligation: 9 μ l
146 nuclease-free water, 2 μ l 10 mM dNTPs, 8 μ l 5x SuperScript IV first
147 strand buffer, 4 μ l 0.1 mM DTT.
- 148 3. Add the reverse transcription master mix to the ligation reaction and
149 mix by pipetting. Add 2 μ l SuperScript IV reverse transcriptase and
150 mix by pipetting. Incubate in a thermocycler with the following pro-
151 tocol: 50 min at 50 °C, 10 min at 70 °C, cool down to 4 °C.
- 152 4. Let the Agencourt RNAClean XP beads come to room temperature
153 during reverse transcription. Carefully resuspend beads before use.
154 Transfer reaction to a 1.5 ml DNA LoBind tube and mix with 72 μ l
155 Agencourt RNAClean XP beads. Incubate 5 min at room temperature
156 on a gentle rotator mixer.
- 157 5. Collect beads on a magnetic stand and remove supernatant. Wash
158 pelleted beads two times (30 sec) with 200 μ l freshly prepared 70 %
159 ethanol. Remove supernatant. Spin sample down and place on magnet
160 again. Remove any residual ethanol.
- 161 6. Resuspend beads in 20 μ l nuclease-free water by gentle flicking and
162 incubate 5 min at room temperature on a gentle rotator mixer. Collect
163 beads on a magnetic stand and transfer 20 μ l eluate in a fresh 1.5 ml
164 DNA LoBind tube.
- 165 7. For ligation of the RMX adapter, add the following to 20 μ l eluate: 8
166 μ l NEBNext Quick Ligation Reaction Buffer, 6 μ l RMX (from SQK-
167 RNA002), 3 μ l nuclease-free water, 3 μ l T4 DNA Ligase. Mix by
168 pipetting and incubate 10 min at room temperature.
- 169 8. Add 40 μ l carefully resuspended Agencourt RNAClean XP beads to
170 the reaction and mix by pipetting. Incubate 5 min at room tempera-
171 ture on a gentle rotator mixer.

- 172 9. Collect beads on a magnetic stand and remove supernatant. Wash
173 pelleted beads two times with 150 μ l wash buffer (WSB, from SQK-
174 RNA002). Resuspend beads by flicking, spin down and return to mag-
175 netic stand. Remove supernatant from pelleted beads.
- 176 10. Resuspend beads in 21 μ l elution buffer (EB, from SQK-RNA002) by
177 gentle flicking and incubate 5 min at room temperature on a gentle
178 rotator mixer. Pellet beads on a magnetic stand and transfer 21 μ l
179 eluate in a fresh 1.5 ml DNA LoBind tube.
- 180 11. Quantify 1 μ l of the library on a Qubit fluorometer with the Qubit
181 dsDNA HS kit according to the manufacturerers protocol. Concentra-
182 tion should be usually in the range of 5 - 10 ng/ μ l.
- 183 12. Insert MinION R9.4.1 Flow cell in the MinION or GridION sequenc-
184 ing device and perform Flow cell check in the MinKNOW software.
185 For successful sequencing of mammalian polyA⁺ RNA at least 1,000
186 available pores are recommended.
- 187 13. Prepare Priming Mix by adding 30 μ l flush tether (FLT, from EXP-
188 FLP002) to a vial of flush buffer (FB, from EXP-FLP002) and mix by
189 pipetting. Open priming port. Remove air bubble from priming port
190 by inserting the tip of a P1000 pipette into the priming port and slowly
191 dialing up, until a small volume of storage buffer enters the pipette
192 tip. Load 800 μ l Priming Mix via the priming port and carefully avoid
193 introduction of air bubbles. Close the priming port and wait for 5 min.
- 194 14. Mix 20 μ l library with 17.5 μ l nuclease-free water and 37.5 μ l RNA run-
195 ning buffer (RRB, from SQK-RNA002) and mix by pipetting. Open
196 the priming port and the sample port. Load 200 μ l Priming Mix via
197 the priming port. Mix library by pipetting just before loading and
198 load dropwise via the sample port. Carefully avoid introduction of air
199 bubbles. Close the sample port and the priming port.
- 200 15. Start sequencing for 48 to 72 h in the MinKNOW software. Choose
201 direct RNA-sequencing kit and high-accuracy basecalling as paramet-
202 ers. We recommend to adjust the output filter to a minimum Q score
203 of 7 (instead of 9).

204 **Preparation of an *in vitro* transcriptome sample**

205 The *in vitro* transcriptome sample is prepared based on a protocol published
206 by Zhang *et al.* ? with some modifications.

- 207 1. Adjust 100 ng polyA⁺ RNA to a total volume of 6 μ l with nuclease-
208 free water. Add 1 μ l each of 10 μ M oligo(dT)-VN RT primer and 10

- 209 mM dNTPs. Mix by pipetting and incubate in a thermocycler: 5 min
210 at 75 °C, 2 min at 42 °C, cool to 4 °C.
- 211 2. Assemble 2.5 µl 4x template switching RT buffer, 0.5 µl 20 µM TSO,
212 1 µl 10x template switching RT enzyme mix and mix by pipetting.
213 Combine with 6 µl RNA and incubate in a thermocycler: 90 min at
214 42 °C, 10 min at 68 °C, cool to 4 °C.
- 215 3. For Second strand synthesis add to First strand synthesis reaction: 50
216 µl Q5 Hot Start High-Fidelity 2X Master Mix, 5 µl RNase H, 2 µl 10
217 µM T7 extension primer, 33 µl nuclease-free water. Mix by pipetting
218 and incubate in a thermocycler: 15 min at 37 °C, 1 min at 95 °C, 10
219 min at 65 °C, cool to 4 °C.
- 220 4. Purify double stranded cDNA with NucleoSpin Gel and PCR Clean-up
221 kit according to the manufacturer's protocol and elute in 20 µl elution
222 buffer. Determine concentration on a Nanodrop spectrophotometer.
223 cDNA may be stored at -20 °C.
- 224 5. Combine 8 µl cDNA for *in vitro* transcription with 2 µl each of ATP,
225 GTP, CTP, UTP, 10x reaction buffer and enzyme mix from the MEGAscript
226 T7 transcription kit. Incubate 3 h at 37 °C.
- 227 6. Digest template DNA by addition of 1 µl Turbo DNase. Mix by pipet-
228 ting and incubate 15 min at 37 °C.
- 229 7. Adjust reaction volume to 100 µl with nuclease-free water and clean up
230 with RNA Clean & Concentrator-25 kit according to the manufactur-
231 ers protocol, using two volumes of adjusted RNA binding buffer (1:1
232 RNA binding buffer : ethanol). Elute RNA in 25 µl nuclease-free wa-
233 ter. Determine RNA concentration on a Nanodrop spectrophotometer.
234 Store at -80 °C.

235 Nanopore read processing

236 Minimap2 and samtools

237 Christoph

238 Use Case 1: Comparison of wildtype and knock-out samples

239 Xpore data

240 Christoph

241 Use Case 2: Comparison of wildtype and IVT samples

242 Christoph

243 **Use Case 3: Comparison of wildtype to simulated IVT sample**

244 **Christoph**

245 **NOTES**

246 **Tips and Tricks**

247 **ACKNOWLEDGMENTS**

248 The authors would like to thank Etienne Boileau, Thiago Britto Borges,
249 Tobias Jakobi for proof-reading and comments. The authors are grateful
250 to Marek Franitza for running the experiments on the 10x platform and to
251 Christian Becker for running ONT sequencing. This work was supported by
252 Informatics for Life funded by the Klaus Tschira Foundation.

FIGURE CAPTIONS




Figure 1: **General outline of RNA modification detection by JACUSA2.** A key feature of our approach is that multiple replicates can be compared as shown on the left. Samples of interests where all modifications are present could be compared with either KO samples where the modification of interest is missing or IVT/simulated samples where all modifications are absent. Read stacks (in blue) are compared head-to-head as shown on the right.




Figure 2: **Motivation of 5mer context for RNA modification mapping.** The nanopore covers 5 consecutive RNA residues. That is why we consider a 5mer context and derive 3 principal features for every position within a given 5 mer (15 features in total, with a central A residue in this example). We evaluate each feature set by previously learned patterns and compute a final score for modification site detection.




Figure 3: **Experimental and computational workflow.** tbd

Software	Version	Description
Minimap2	https://github.com/lh3/minimap2 v2.22 or later	https://lh3.github.io/minimap2/
samtools	https://github.com/samtools/samtools v1.12 or later	http://samtools.github.io/
JAVA	openjdk 11.0.12 2021-07-20 - JAVA 11 or later	OpenJDK Runtime Environment
R	https://www.r-project.org/ version 3.5.1 or later	The R Project for Statistical Computing
PERL	https://www.perl.org/ version 5.28.1 or later	Perl is a highly capable, feature-rich programming language
BASH, sed, awk	should be part of your Linux distribution	Misc.
bedtools	https://github.com/arq5x/bedtools2 version 2.29.2 or later	Perl is a highly capable, feature-rich programming language
NanoSim	https://github.com/bcgsc/NanoSim version 3.0.2 or later (optional)	NanoSim is a fast and scalable read simulator that captures the technology-specific features of ONT data

Table 1: **Software dependencies** blubba

254 TABLE CAPTIONS

255 TABLES

R Pack- ages	Version	Description
ggplot2	https://cran.r-project.org/web/packages/ggplot2/index.html - ggplot2_3.3.0 or later	ggplot2 is a system for declaratively creating graphics, based on The Grammar of Graphics.
NMF	https://cran.r-project.org/web/packages/NMF/index.html - NMF_0.22.0 or later	Provides a framework to perform Non-negative Matrix Factorization (NMF).

Table 2: **R Package dependencies** blubba